



Visual Information Processing
Department of Computing, Imperial College

Imperial College
London




A Stereoscopic Fibroscope for Camera Motion and 3D Depth Recovery during Minimally Invasive Surgery

**David Noonan, Peter Mountney, Daniel Elson,
Ara Darzi and Guang-Zhong Yang**

Institute of Biomedical Engineering & Dept. of Biosurgery & Surgical Technology
Imperial College London

RobotsSLAM in Minimally Invasive Surgery

- Use of robotic surgery systems increases
- Scope exists to enable ever more complex procedures
- Interventions along curved instrument pathways
- Stereo fibre image guide based system is ideally placed for integration with flexible instrumentation where miniaturisation is required






Outline


- System design of a Stereoscopic Fibroscope
- SLAM algorithm design
- Validation steps performed
- Results
- Additional capabilities

Fibre Image Guides

- Type: Sumitomo IGN-05/10 coherent fibre image guide
 - Pixels: 10,000
 - Length: 1.5mm
 - Diameter: 590 μ m
 - Min. bending radius: 25mm
- GRIN lens cemented to each tip
 - Image area: 35x35mm²
 - Working distance 20mm



[Cross Sectional View of Image Guide Fiber]



Number of picture elements	10000
Jacketing diameter (um)	500
Picture elements area diameter (um)	450
Coating diameter - primary (um)	590
Coating diameter – secondary (um)	---
Circularity	>/=0.93
Core material	GeO2 containing silica
Cladding material	F containing silica Silicone
Coating material	Silicone
Numerical aperture	0.35
Lattice defect (%)	</=0.1
Allowable bending radius	25
Allowable max temp (°C)	150



Stereoscopic Fibroscope

1. Passive 3-axis joint to facilitate free-hand motion
2. Rigid body to mount optical tracking markers
3. Protective tubing for fibre image guides
4. 2 * 10,000 flexible fibre image guide
5. Grubscrew to adjust camera vergence
6. Tubing path to image acquisition system






Image Acquisition System




- Acquisition via:
 - CCD Camera: UEye, UI-2250-C/CM
 - Both images focussed onto a single camera
 - X10 achromatic microscope objective
 - 100mm focal length lens
 - XYZ positioning stage
 - Custom fibre mount

Full System





Full System (2)






Calibration





- Left camera image co-ordinate system defined manually
 - Due to arbitrary rotation of image around camera z-axis
- Right camera image co-aligned with that of the left image
- Stereo camera calibration performed
 - Camera intrinsic & extrinsic parameters
 - Manual corner detection due to poor image resolution
- Hand-Eye calibration
 - Rigid body origin to left camera centre



Preprocessing




- Honeycomb artifact removal
 - Honeycomb fibre structure visible
 - Bandpass filter in frequency space and Gaussian smoothing







SLAM Algorithm Design

- **Initialise: Feature detection**
- **Triangulate: 3D map feature**
- **Estimate Camera Motion**
- **Match Features: Active Search**
- **Update Camera and Map**







Validation Setup



Results (1)







Camera Motion



Results (3)

Surface Reconstruction

- Overall reconstruction error:
 - X – 2mm
 - Y – 1.3mm
 - Z – 2.9mm
- Z error larger due to small baseline




CT Ground Truth

SLAM reconstruction





Results (2)

Camera Motion



Conclusion

- Demonstrated feasibility of integrating twin flexible fibre guides in stereo configuration
- Described the challenges
- Successfully employed within a SLAM algorithm
 - Motion recovery
 - 3D Depth recovery
- Potential application
 - Within a catheter which utilizes the stereo vision for targeting and depth recovery for accurate focussed energy delivery






Acknowledgments

wellcometrust

EPSRC Engineering and Physical Sciences
Research Council



The Hamlyn Centre for Robotic Surgery

Supported by Helen Hamlyn Trust

Technology Strategy Board

Driving Innovation

Robotic Assisted Surgery

- Use of robotic surgery systems increases
- ≈ 900 units installed worldwide
- Majority of cases for radical prostatectomy
- Scope exists to enable ever more complex procedures

